

# Phantom

SAFETY, RELIABILITY  
& COMFORT



Protection from the harshest of weather conditions with a premium impact resistant, clear vision canopy. Multiple tiller settings and a plush high-back captain's chair ensure a comfortable driving position. The Phantom comes standard with a premium lighting package coupled with incredible storage capacity in front and rear baskets. High profile tyres synced with sensational all round suspension deliver superior driver comfort. Australia's best value for money canopy scooter!



**TopGun**  
**Mobility**  
WHEELS FOR LIFE

## Phantom Features

- Sturdy canopy
- Front protective windscreen
- Seat belt
- Deluxe orthopedic reclining seat
- Angle adjustment fold up arm rests
- Automatic braking system
- Intelligent LCD digital display
- LED headlight
- Direct drive sealed transmission/gearbox
- Multi position tiller angle
- Free-wheeling device
- USB charging
- Rear basket
- Front storage box with 2 drink holders
- Power save mode
- 5 year limited warranty

## Specifications

Motor Power	950W
Max. Speed	10km/h
Operation Range	Up to 50km
Max. Climbing Grade	10° (depending on weight of user)
Turning Radius	185cm
Ground Clearance	15.5cm (+/- 2.5cm)
Max. Loading Weight	180kg
Overall Dimensions	130 x 83 x 168cm
Wheel Size - Front	13" pneumatic
Wheel Size - Rear	13" pneumatic
Weight	96kg
Battery	24V 70Ah x 2 (lead-acid)
Controller	PG 140A
Brake System	Electromagnetic
Free Wheel Mode	Yes
Suspension	Yes
Rear Mirrors	Left & Right
Charger Size	24V 8A



**NOTE: All specifications subject to change without notice.**

- The information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery Condition and tyre pressure.