



LIGHTWEIGHT ALUMINIUM WHEELCHAIR

Product features

- Constructed from sturdy aluminium, this wheelchair is durable and robust, yet lightweight for transporting
- Nylon seat and back offer comfortable support and are easy to clean
- Chair folds up easily for storage and transportation, further enhanced by half-folding backrest, quick-release footrests and rear wheels (self-propelling models only)
- Full-length armrests with upholstered arm pads increase user comfort and are easy to wipe clean
- Footrest height can be adjusted tool-free into one of 5 positions to optimise user comfort and support
- Footrests can also be swung out of the way to facilitate transfers on and off chair
- Rear stepper tubes allow carers/ attendants to safely climb curbs
- Puncture-proof front wheels ensure low maintenance and longevity
- Supplied with lap belt for increased user safety



Self-Propel Model

- Aluminium 61cm (24") rear wheels with hand rim and grey pneumatic tyres allow a smooth independent ride for active users whilst making chair lightweight for transporting
- Total weight: only 15kg (33lb)

Transit Model

- Cable brakes on rear handles allow immediate access for attendant to reduce speed or lock brakes for parking
- Total weight: only 12kg (27lb) for easy transportation

“The Lightweight Aluminium Wheelchair offers a value-for-money introduction to the Aluminium Wheelchair range without compromising on quality.”

Code	Width (Open)	Width (Closed)	Length	Height	Seat Width	Seat Length	Seat Height	Weight	Weight Capacity
LAWC001AU	61cm (24")	24cm (9")	106cm (42")	90cm (35")	46cm (18")	40cm (16")	48cm (19")	15kg (33lb)	115kg (18st)
LAWC002AU	55cm (22")	29cm (11")	96cm (38")	91cm (36")	46cm (18")	40cm (16")	48cm (19")	12kg (27lb)	115kg (18st)

“This Ultra Lightweight Wheelchair combines style with function and practicality.”

ULTRA LIGHTWEIGHT ALUMINIUM WHEELCHAIR



Product features

- Ultra-lightweight aluminium frame with single cross brace is strong and stable, yet very lightweight for transporting and easily folded for storage
- Available in 18" and 20" seat widths
- Half-folding backrest, quick-release footrests and rear wheels (self-propelling model only) further facilitates storage and transportation
- Grey two-tone upholstery offers comfortable support, is easy to clean and along with silver frame has an attractive appearance
- Soft-touch fabric side panels improve comfort
- Features handy rear pocket for storing items
- Footrests easily adjust into one of five positions to suit user support requirements and swing out of the way to facilitate transfers
- Desk-style armrests allow easy access to work surfaces
- Two stepper tubes enable attendant to safely raise front wheels for mounting curbs
- Two solid Polyurethane (PU) front tyres, 20cm (8") in diameter are durable, low maintenance and assist with manoeuvrability
- Long-reach brakes on rear tyres are easily activated by user or attendant for stable parking
- Spoked 61cm (24") rear wheels with solid grey Polyurethane (PU) tyres are hard-wearing and low maintenance and together with lightweight aluminium rims and hand rims increase manoeuvrability (Self Propel Model)
- Cranked push handles allow for compact folding (Transit Model)
- Cable brakes on rear handles provide attendant with immediate access to reduce speed or to lock for stable parking (Transit Model)
- 32cm (12.5") solid rubber rear wheels assist with manoeuvrability and are durable and low maintenance (Transit Model)



Code	Width (Open)	Width (Closed)	Length	Height	Seat Width	Seat Length	Seat Height	Weight	Weight Capacity
LAWC007AAU	62cm (24")	32cm (13")	107cm (42")	93cm (37")	42cm (17")	39cm (15.5")	52cm (20.5")	13kg (29lb)	115kg (18st)
LAWC008AAU	62cm (24")	28cm (11")	95cm (37")	97cm (38")	42cm (17")	39cm (15.5")	52cm (20.5")	12kg (26lb)	115kg (18st)
LAWC011AAU	71cm (28")	32cm (13")	109cm (43")	93cm (37")	51cm (20")	39cm (15.5")	52cm (20.5")	13kg (29lb)	125kg (20st)
LAWC012AAU	66cm (26")	28cm (11")	95cm (37")	97cm (38")	51cm (20")	39cm (15.5")	52cm (20.5")	12kg (26lb)	125kg (20st)